# LEPL1201 - Laboratoire 1

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### 3 Instructions (from the field mapper kit manual)

**Important :** the silver conductive ink reaches its maximum conductivity after 20 minutes drying time. For optimal results plan the timetable for conducting the experiments and correlate drawing the conductive ink paths accordingly.

1. Plan and sketch the layout (size, shape and relative spacing) of the charged paths to be studied on a piece of scratch paper. These paths can be any two dimensional shape, such as straight or curved lines, circles, dots, squares, etc. Since the charged paths will actually be conductive ink electrodes, they will be referred to as electrodes.
2. Draw the electrodes on the black paper. This is the most difficult and crucial part of the experiment. Follow these steps carefully.
  - Place the conductive paper, printed side up, on a smooth hard surface. Do not attempt to draw the electrodes while the paper is on the corkboard.
  - Shake the ink bottle vigorously for 10-20 seconds to disperse any particle matter suspended in the ink.
  - Draw the electrodes on the black conductive paper. If the line becomes thin or spotty, draw over it again. A solid line is essential for good measurements. The line will be air

dry in 3-5 minutes at room temperature. However, the medium won't reach maximum conductivity until after 20 minutes drying time.

3. Mount the conductive paper on the corkboard using one of the metal push pins in each corner.
4. Connect the electrodes to a battery, DC power supply, or any other potential source in the 5 to 20 VDC range using the supplied connecting wires (see figure 3a). The potential source should be capable of supplying 25 mA. (if possible, the potential should be equal to the full scale reading of the electronic voltmeter used in the experiment)
  - Place the terminal of a connecting wire over the electrode, then stick a metal push pin through its terminal and the electrode into the corkboard. Make certain that the pin holds the terminal firmly to the electrode (see figure 3b).

Note : check that the surface of the terminal which touches the electrode is clean. A dirty path may result in a bad contact.

Connect the other end of the wire to the battery.

5. To check the electrode proper conductivity connect one voltmeter lead near the push pin on an electrode. Touch the voltmeter's second lead to other points on the same electrode. If the electrode has been properly drawn, the maximum potential between any two points on the same electrode will not exceed 1% of the potential applied between the two electrodes.

Note : this test can only be made if the potential source is connected across the two electrodes.

If the voltage across the same electrode is greater than 1% of the voltage applied between the two electrodes, then remove the paper from the corkboard and draw over the electrodes a second time with the conductive ink.

6. Equipotentials are plotted by connecting one lead of the voltmeter (the ground) to one of the electrode push pins. This electrode now becomes the reference. The other voltmeter lead (the probe) is used to measure the potential at any point on the paper simply by touching the probe to the paper at that point.

To map an equipotential, move the probe until the desired potential is indicated on the voltmeter. Mark the paper at this point with a soft lead or light-colored lead pencil. Continue to move the probe, but only in a direction which maintains the voltmeter at the same reading. Continue to mark these points. Connecting the points produces an equipotential line.